

Omicron Variant Performance Study Report

COVID-19 Antigen Detection Kit

Sponsor:

New Gene (Hangzhou) Bioengineering Co., Ltd.

Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street

Binjiang District

Hangzhou City, Zhejiang Province

P. R. China

<https://www.new-gene.net/>

Content

SPONSOR - INVESTIGATION -STUDY COORDINATION.....	3
SPONSOR	3
INVESTIGATION	3
STUDY COORDINATION DEPARTMENT	3
1. BACKGROUND.....	4
2. PURPOSE	5
3. PRODUCT INFORMATION:	5
4. PRINCIPLE.....	7
5. STUDY DETAILS	8
6. EVALUATION SITE DETAILS	9
7. TEST RESULT.....	10
8. RESULT ANALYSIS.....	11
9. CONCLUSION.....	12
10. SIGNATURES&APPROVAL	12
11. REFERENCE.....	13
12. ANNEX I - LABORATORY VERIFICATION DATA LIST	14

SPONSOR - INVESTIGATION - STUDY COORDINATION

SPONSOR:

New Gene (Hangzhou) Bioengineering Co., Ltd.

Room 1606, Floor 16, Building 5, 688 Bin'an Road, Changhe Street

Binjiang District

Hangzhou City, Zhejiang Province

P. R. China

(+86) 0571-56515020

marketing@new-gene.com

<https://www.new-gene.net/>

INVESTIGATION:

iLab Medical AB

Järnbrotts Prästväg 2,

421 47 Västra Frölunda

Sweden



STUDY COORDINATION DEPARTMENT:

VIKBAY AB

Maj På Malös Gata 21

417 67, Göteborg, Sweden

Tel.: (+46)7347-87389

e-mail: adaline.lee@vikbay.se



Homepage: www.vikbay.com

© 2022 VIKBAY. All rights reserved

Evaluation of COVID-19 Antigen Detection Kit For Omicron Variant

New Gene (Hangzhou) Bioengineering Co., Ltd.

1. BACKGROUND

The COVID-19 pandemic, caused by the severe acute respiratory syndrome coronavirus 2 virus (SARS-CoV-2), has been ravaging the world for more than two years. Many efforts have been made to curb the spread of SARS-CoV-2 in societies. Rapid and accurate detection of SARS-CoV-2 is crucial to limit virus transmission. Antigen immunochromatographic tests have been widely used in point of care settings due to it is a fast, cheap, easy-to-use method [1]. However, new variants of SARS-CoV-2, mainly variants of concern (VOC), with increased transmissibility, have brought additional challenges to the current prevention and control measures of SARS-CoV-2. Recently, omicron variant of SARS-CoV-2 has rapidly replaced delta variant in a short period in many European countries [2], the daily omicron infection case is more than ten times the peak of delta wave [3]. Numerous mutations have been identified in the receptor-binding domain of spike glycoprotein of the omicron variant [4], which raised the concerns about the performance of current rapid antigen kits that are solely targeting the spike protein of omicron variant may be impaired (thus not combined with the nucleocapsid protein) as well as the viral load measured at different time points and at different sites (e.g. throat and nose) after Omicron infection. Therefore, an evaluation of the performance of rapid antigen tests for the omicron variant is needed [5].




Figure 1. Omicron genome of SARS-CoV-2 with accumulations of mutations in S-Gene [6].

2. PURPOSE

The objective of this study was to evaluate the performance of the New Gene COVID-19 Antigen Detection Kit for the detection of the omicron variant of SARS-CoV-2.

3. PRODUCT INFORMATION:

3.1 Detection System

Manufacturer Name	New Gene (Hangzhou) Bioengineering Co., Ltd.
Test Name	COVID-19 Antigen Detection Kit
Pack Size(s)	25 tests/kit
Contents of Kit	Aluminum Foil Pouch, Test Card, Desiccant, Instructions for User, Sampling Swab, Sample Extraction Tube & Tube Cap
Lot No.	20211120-01
Prod. Date	2021-11-20
Exp. Date	2023-11-19
Product Storage (temperature range)	2-30°C
Shelf-life (months)	24 months
Manufacturing Site (country)	Hangzhou, China

3.2. Reference PCR Method

Reference Method	TaqPath™ COVID-19 CE-IVD RT-PCR Kit		
Pack Size(s)	1,000 reactions	Shelf-life (months)	12 months
Cat No.	A48067	Manufacturer	Thermo Fisher Scientific Inc.
Storage	-30°C to -10°C; TaqPath™ COVID-19 Control, ≤ -70°C keep away from light		
Product Name	Product No.	Lot No.	Expiration Date
TaqPath™ COVID-19 Control	A48105	248106	16-Dec-22
TaqPath™ COVID-19 Control Dilution Buffer	A48099	2201105	14-Jan-23
TaqPath™ 1-Step Multiplex Master Mix	A48111	01206471	20-Dec-22
TaqPath™ COVID-19 RT-PCR CE-IVD KIT	A48012	2201238	20-Dec-22

3.3 Variant of Concern Reagents Used in the Evaluation

Reference Method	Omicron-specific PCR test for new Variant of Concern		
Pack Size(s)	96-10.000 reactions	Catalog Number	031522
Expiration Date	2023-01-30	Storage	-30°C to -10°C
Lot No.	C52-32-32	Shelf-life (months)	12 months
Contents of kit	GenXPro Covid-2 (Delta, Omicron) detection kit		
Manufacturer	GenXPro GmbH, Altenhöferallee 3, 60438 Frankfurt am Main, Germany		

4. PRINCIPLE

The New Gene COVID-19 Antigen Detection Kit is an immunochromatographic membrane assay that uses highly sensitive monoclonal antibodies to detect the nucleocapsid protein of SARS-CoV-2 in nasal swab. The device is composed of four parts: sample pad, reagent pad, reaction membrane, and absorbing pad. The reagent pad contains the colloidal-gold conjugated with the monoclonal antibody against the nucleocapsid protein of SARS-CoV-2. The reaction membrane contains the secondary antibodies for nucleocapsid protein of SARS-CoV-2. When the sample is added into the sample well, conjugates absorbed in the reagent pad are dissolved and migrate along with the sample. If SARS-CoV-2 antigen is present in the sample, the complex of the anti-SARS-CoV-2 conjugate and the virus will be captured by the specific anti-SARS-CoV-2 monoclonal antibodies coated on the test line region (T). The absence of the T line suggests a negative result. To serve as a procedural control, a red line will always appear in the control line region (C) indicating that proper volume of sample has been added and membrane wicking effect has occurred.

5. STUDY DETAILS

Study Design:	<p>Prospective performance evaluation studies were sampled at multiple independent sites in Sweden. The purpose of this study was to establish the performance of New Gene COVID-19 Antigen Detection Kit for the omicron variant. Presence of symptoms, duration of symptoms were collected for all enrolled samples.</p> <p>The nasal swab sample was tested with the New Gene COVID-19 Antigen Detection Kit according to the manufacturer's instructions in iLab Medical Laboratory. The same NPS sample was tested by reference RT-PCR method in same laboratory. The result of the New Gene COVID-19 Antigen kit was compared with the result of the RT-PCR assay to validate the performance of the New Gene COVID-19 Antigen Detection Kit for the detection of the omicron variant of SARS-CoV-2.</p>
fFEvaluation Period:	20220121 to 20220206
Sample Type:	Nasal swabs
Inclusion Criteria:	<ol style="list-style-type: none"> 1. Individuals were ≥ 18 years old; 2. Samples from individuals who were symptomatic and suspected of COVID-19, as well as asymptomatic individuals; 3. The RT-PCR test on nasopharyngeal swabs was positive for SARS-CoV-2. ; 4. At all ages, male or female (try to cover all age groups).
Exclusion Criteria:	<ol style="list-style-type: none"> 1. Sample collection, shipment and storage do not meet relevant requirements; 2. Abnormal experimental results caused by operation errors and/or kit failure.
Sample Size:	A total of 150 nasopharyngeal swabs from healthcare settings were registered in this study. Among them, 41 samples were tested SARS-CoV-2 positive by RT-PCR.

Reference Method:	RT-PCR assay was used as the reference method. The RNA was extracted from nasopharyngeal swabs using Total RNA Purification 96-Well Kit (Norgen Biotek, Canada), and the TaqPath™ COVID-19 CE-IVD RT-PCR Kit (Thermo Fisher Scientific, USA) was used for RT-PCR test. The run was on Roche® 480 Light Cycle II system (Roche Diagnostics, Switzerland).
Laboratory Verification Performance	<ol style="list-style-type: none"> 1. Sensitivity was calculated as the proportion of samples tested true positive by New Gene COVID-19 Antigen Detection Kit among all positive samples tested by the RT-PCR method. The sensitivity is reported as a percentage. 2. Specificity was not available in this study. 3. The 95% confidence intervals (CI) were calculated to assess the level of uncertainty introduced by sample size.

6. EVALUATION SITE DETAILS

Country	Sweden
Location of Laboratory Verification Site(s) (city, town)	Gothenburg
Health Care Level of Site(s)	Community Testing Clinic
Sample Collection	iLab Medical AB
Research Unit	iLab Medical AB
Address	Järnbrotts Prästväg 2, 421 47 Västra Frölunda
Laboratory Validation (Variant of Concern Test)	iLab Medical

7. TEST RESULT

7.1 Study Cohort

Country	Sweden
Total N (valid PCR results)	41
Age [mean (min-max), N]	31.32 (19-50), 41
Gender [%F, (n/N)]	41.46% (17/41)
Positivity [%, (n/N)]	97.56% (40/41)
Symptoms Present¹ [%Yes, (n/N)]	21.95% (9/41)
Days From Symptom Onset¹	
Days < 0-3 (n, %)	7, 77.8%
Days 4-7 (n, %)	0%
Days 8+ (n, %)	2, 22.2 %
PCR Ct [median; N]	23; 41
Ct ≥ 30 (n, %)	5, 12%
25 < Ct < 30 (n, %)	8, 20%
Ct ≤ 25 (n, %)	28, 68%

¹Note: data on the onset of symptoms are only available for individuals who tested positive for PCR..

8. RESULT ANALYSIS

8.1. Laboratory Verification Performance against the Reference Method

The performance of New Gene COVID-19 Antigen Detection Kit was established with 41 nasal swabs collected from symptomatic individuals who were suspected of COVID-19 as well as asymptomatic individuals..

New Gene COVID-19 Antigen Detection Kit	Reference RT-PCR Method		
	+	-	Total
+	40	0	40
-	1	0	1
Total	41	0	41
Sensitivity = 40/41 = 97.6% (95% CI: 87.1% ~ 99.9%)			
Specificity Not available			

8.2 Laboratory Verification Performance against the Reference Method – by Cycle Threshold Counts

The performance of the New Gene COVID-19 Antigen Detection Kit with positive results stratified by the reference assay cycle threshold (Ct) counts was assessed. As presented in the table below, the positive agreement of the New Gene COVID-19 Antigen Detection Kit is higher with samples of a Ct count < 30.

New Gene COVID-19 Antigen Detection Kit	Reference RT-PCR Method		
	Pos. (Ct ≤ 25)	Pos. (25 < Ct < 30)	Pos. (Ct ≥ 30)
Positive	28	8	4
Negative	0	0	1
Total	28	8	5
Positive Agreement (95% CI)	100% (87.9% ~ 100%)	100% (67.6% ~ 100%)	80% (37.6% ~ 96.4%)

9. CONCLUSION

The evaluations showed that New Gene COVID-19 Antigen Detection Kit has a good agreement with reference RT-PCR Method for omicron variant. The laboratory verification sensitivity of New Gene COVID-19 Antigen Detection Kit is 97.6%, and it is 100% sensitive for samples with Ct value < 30.

All the New Gene COVID-19 Antigen Detection Kit tested in this work perform well. The tests show great sensitivity. The laboratory evaluation might differ from the performance on the field due to the limited sample range and close monitoring is needed by professionals. From the results of 41 tested samples, it presents that the New Gene COVID-19 Antigen Detection Kit developed by New Gene (Hangzhou) Bioengineering Co., Ltd. meets the requirement of intended use.

10. SIGNATURES&APPROVAL

Evaluation site: Göteborg, Västra Götalands län, Sweden

Principal evaluator: Shokoofeh Naghdipour, Dr. Huaqing Li

Reported by:

 _____

Date:

2022-02-10

Signature:

 _____

Stamp:



11. REFERENCE

1. Kubina R, Dziedzic A: **Molecular and Serological Tests for COVID-19 a Comparative Review of SARS-CoV-2 Coronavirus Laboratory and Point-of-Care Diagnostics.** *Diagnostics* 2020, **10**(6).
2. **Weekly epidemiological update: Omicron variant of concern (VOC) – week 2 (data as of 20 January 2022) EU/EEA** [<https://www.ecdc.europa.eu/en/news-events/weekly-epidemiological-update-omicron-variant-concern-voc-week-2-data-20-january-2022>]
3. Murray CJ: **COVID-19 will continue but the end of the pandemic is near.** *The Lancet* 2022.
4. Saxena SK, Kumar S, Ansari S, Paweska JT, Maurya VK, Tripathi AK, Abdel-Moneim AS: **Characterization of the novel SARS-CoV-2 Omicron (B.1.1.529) variant of concern and its global perspective.** *J Med Virol* 2021.
5. Joint Research Centre (2022, January 26). **Additional information for manufacturers.**
Retrieved from <https://covid-19-diagnostics.jrc.ec.europa.eu/manufacturers-info>
6. GenXPro GmbH. (2022). **Omicron-specific PCR test for new Variant of Concern.** Retrieved from <https://genxpro.net/sars-cov-2-omicron-pcr-test/>

12. ANNEX I - LABORATORY VERIFICATION DATA LIST

No.	Gender	Age	Sample Collection Date	Symptom Days	Symptoms when specimen collected: A: please describe the symptom (like Cough, Fever, Dizzy...) B: Asymptomatic	Result									
						Antigen Rapid Test				PCR Result					
						Antigen Lot Number	Specimen Type	New Gene Testing Date	Antigen test result	PCR Confirmation Date	Specimen Type	ORF1ab CT Value	SARS-CoV-2 variant type	PCR operation date	Result
Omi01	M	20	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	21.06	Omicron	20220204	Positive
Omi02	F	33	20220204	0-3	Cough, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	20.59	Omicron	20220204	Positive
Omi03	M	29	20220204	8+	Cough, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	21.67	Omicron	20220204	Positive
Omi04	M	46	20220204	0-3	Cough, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	21.21	Omicron	20220204	Positive
Omi05	F	42	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	20.95	Omicron	20220204	Positive
Omi06	F	19	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	20.62	Omicron	20220204	Positive
Omi07	M	36	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.17	Omicron	20220204	Positive
Omi08	M	50	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	20.00	Omicron	20220204	Positive
Omi09	M	44	20220204	0-3	Cough, Dizzy	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	21.29	Omicron	20220204	Positive
Omi10	M	32	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	20.93	Omicron	20220204	Positive
Omi11	F	24	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.64	Omicron	20220204	Positive

Omi12	F	26	20220204	8+	Dizzy, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.78	Omicron	20220204	Positive
Omi13	M	43	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.15	Omicron	20220204	Positive
Omi14	M	22	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	23.00	Omicron	20220204	Positive
Omi15	M	41	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.21	Omicron	20220204	Positive
Omi16	M	27	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.35	Omicron	20220204	Positive
Omi17	M	28	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.51	Omicron	20220204	Positive
Omi18	M	39	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	21.49	Omicron	20220204	Positive
Omi19	M	28	20220204	0-3	Cough, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.14	Omicron	20220204	Positive
Omi20	M	29	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	24.03	Omicron	20220204	Positive
Omi21	M	32	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.56	Omicron	20220204	Positive
Omi22	M	29	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	23.44	Omicron	20220204	Positive
Omi23	M	28	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	23.84	Omicron	20220204	Positive
Omi24	F	28	20220204	0-3	Cough, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	24.49	Omicron	20220204	Positive
Omi25	F	27	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	22.54	Omicron	20220204	Positive
Omi26	M	29	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	23.89	Omicron	20220204	Positive
Omi27	F	28	20220204	0-3	Dizzfy, Fever	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	25.03	Omicron	20220204	Positive

Omi28	F	27	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	23.88	Omicron	20220204	Positive
Omi29	F	35	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	23.49	Omicron	20220204	Positive
Omi30	F	34	20220204	NA	B	20211120-01	Nasal	20220204	Positive	20220204	Nasopharyngeal	25.19	Omicron	20220204	Positive
Omi31	M	22	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	26.26	Omicron	20220204	Positive
Omi32	F	34	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	26.64	Omicron	20220204	Positive
Omi33	F	33	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	25.73	Omicron	20220204	Positive
Omi34	M	39	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	25.83	Omicron	20220204	Positive
Omi35	F	29	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	28.96	Omicron	20220204	Positive
Omi36	M	37	20220204	0-3	Dizzy, Fever	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	27.79	Omicron	20220204	Positive
Omi37	F	33	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	30.2	Omicron	20220204	Positive
Omi38	F	21	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	30.14	Omicron	20220204	Positive
Omi39	F	23	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	30.71	Omicron	20220204	Positive
Omi40	F	24	20220204	NA	B	20211120-01	Nasal	20220206	Positive	20220204	Nasopharyngeal	31.08	Omicron	20220204	Positive
Omi41	M	34	20220204	NA	B	20211120-01	Nasal	20220206	Negative	20220204	Nasopharyngeal	32.66	Omicron	20220204	Positive